

GEROtherm[®] weight system for geothermal probes

for inserting the geothermal probes, materials: cast iron, protection of the geothermal probe and compensation of the uplift force of the geothermal probe, possibility of attaching additional weights. The GEROtherm[®] geothermal probes are mounted without tools using a secure connection with hinge and spring cotter pins.

15 kg initial weight

The weight tip is wedge-shaped, with a socket to connect extension weights, including water guide grooves. Design for tool-free mounting. Length approx. 510 mm Diameter approx. 94 mm Art. No.: 104642 Quantity: Units

19 kg initial weight

The weight tip is wedge-shaped, with a socket to connect extension weights, including water guide grooves. Design for tool-free mounting. Length approx. 615 mm Diameter approx. 94 mm Art. No.: 104643 Quantity: Units

Extension weight 20 kg

The weight tip is wedge-shaped, with a socket to connect additional extension weights, including water guide grooves. Design for tool-free mounting. Length approx. 585 mm Diameter approx. 88 mm Art. No.: 104644 Quantity: Units